



# VACCINATION KITS

- ▶ Our vaccination kits are designed for single use only to avoid cross-infection.
- ▶ Content all necessary supplies for a smart, fast and successful vaccination.
- ▶ Ideal for pandemic preparedness or flu clinic, mobile clinic, off-site vaccination clinic.



Code: VCK01

QTTY.	Description	Instruction For Use
1 ea	Sterile disposable syringe 1-5ML	For withdrawing vaccine from the vial & injecting into recipient
1 pair	Examination gloves	Worn by the health care professional administering the vaccination
1 ea	Drawing up needle 18-21 gauge	For withdrawing the vaccine from the vial
1 ea	Plastic tray	Vessel to contain equipment per recipient.
1 ea	Injecting needle 22-27 gauge	For injecting the recipient. Gauge covers both children & adults, lengths vary.
1 ea	Cotton ball	Placing over the injection site post injection
1 ea	Round band-aid	Placing over the injection site on departure from vaccination clinic
2 ea	Alcohol swab	Wiping the top of the vial + Wiping the injection site (arm) of the recipient
Customization types are possible.		