

# Renal therapy




## Fistula needle



We offer innovative, high performance products and services for use in dialysis. The core of the Renal Therapies spectrum provides a basis which targets improved clinical outcomes for the patient, while minimizing costs and risks for the healthcare providers.


### Features

- ▶ Fixed and rotary wings are available.
- ▶ Japanese made sharp needle proves to be high quality.
- ▶ The large oval orifice on the needle back.
- ▶ Unique design and smooth finish.
- ▶ Super silicon coating of inner and outer needle wall to ease of penetration and minimize the pain.
- ▶ The spiral connector fits perfectly with dialysis line.



## Safety fistula needle

It is used for kidney disease patients to purify the blood, the safety mechanism help to eliminate accidental needle stick injuries.



Code	Gauge	Size, length and color	Description
Non-safety/Safety			
AVFR15/AVFRS15	15G×1"	1.8×25mm Blue	Rotary wing with back eye
AVFR16/AVFRS16	16G×1"	1.6×25mm Green	Rotary wing with back eye
AVFR17/AVFRS17	17G×1"	1.4×25mm Orange	Rotary wing with back eye
AVFF15/AVFFS15	15G×1"	1.8×25mm Blue	Fixed wing with back eye
AVFF16/AVFFS16	16G×1"	1.6×25mm Green	Fixed wing with back eye
AVFF17/AVFFS17	17G×1"	1.4×25mm Orange	Fixed wing with back eye

# Renal therapy



## Blood line



### Features

- ▶ More than 500 types can fit with Fresenius, Cambro, B Braun, Nipro, Baxter hemodialysis machine.
- ▶ Luer type is tight with the adapter of fistula needle.
- ▶ Different clamps fit for tubes of different diameter.
- ▶ Fixed connector with the wing is easy to join the adapter of hemodialysis machine.
- ▶ Customization model are possible upon the requirements.



Code: BL

## Hemodialysis care kit



HCK

Qty	Description	Qty	Description
	Before compartment		After compartment
1ea	Surgical drape 45×45cm	1ea	Surgical drape 45×45cm
1pr	Examination glove	1pr	Examination gloves
5ea	Gauze ball plain type 4"×4"	2ea	Gauze ball plain type 4"×4"
1ea	Alcohol swab	1ea	Alcohol swab
2ea	Adhesive non-woven dressing with open cut in "U" shape	4ea	Gauze swab 7×7cm
4ea	Adhesive non-woven dressing for stabilizing the A-V lines	2ea	Surgical dressing in size 5×8cm
		2ea	Plastic band 40cm

# CATALOGUE CATEGORIES



## Infusion therapy



Page 1-9

- Wide range of a safety way for injection of vein for medication infusion.
- Safety scalp vein set
  - Scalp vein set
  - Infusion set
  - IV burette
  - IV cannula
  - Safety IV cannula
  - IV cannula dressing
  - Sterile IV start kit
  - Extension tube
  - Three way stopcock

## Syringe and needle



Page 10-21

A complete range of syringes for injection or harvesting of biological material with safety protection.

- Disposable syringe with needle
- Disposable syringe with safety needle
- Safety syringe
- Safety needle
- Disposable hypodermic needle
- Auto-Disable syringe for fixed-dose immunization
- Pre-filled flush syringe
- Dental needle
- Irrigation needle
- Suture removal Set
- Veterinary needle
- Hard pack syringe

## Hazardous drug safety



Page 22-28

Various drug preparation needles for safety drug admixture and disposal of used sharps.

- Aspiration needle with safety shield
- Blunt needle with/without filter
- Curved needle
- Side hole needle
- Needle free system
- Safety box for the disposal of used sharps
- Sharp container
- Endoscope transport pad
- Ultimate sharps safety stations (needle counter)
- Safety medical products for NEEDLE STICK PREVENTION

## Safety medical products



Page 28

Wide range of safety medical products to protect health care workers against sharps injuries.

- Safety scalp vein set
- Safety IV cannula
- Safety needle
- Safety syringe
- Auto-Disable syringe
- Safety blood collection needle
- Safety blood lancet
- Safety insulin pen needle
- Safety insulin syringe
- Safety fistula needle
- Safety drug preparation needle
- Ultimate sharps safety stations
- Needle free system
- Safety box for the disposal of used sharps
- Sharp container

## Nutrition therapy



Page 29-32

Various components for nutrition therapy which is compliant with the ISO80369-3 to help reduce the risk of misconnections.

- ENFit@syringe
- The ENFit@ connector
- ENFit@ Feeding set
- Oral syringe
- Irrigation syringe

## Specimen collection



Page 33-36

Our safety blood collection system help eliminate accidental needle stick injuries when collection blood sample.

- Blood collection needle
- Safety blood collection needle
- Safety lancet
- Blood collection tube

## Anesthesia & pain management



Page 37-39

Wide range of products for subdural, lower thorax and lumbar spinal puncture.

- Spinal needle
- Anesthesia kit
- High pressure manifold
- NRFit™ syringe

## Transfusion therapy



Page 40-42

Complete range of safety way to collect, preserve, filtration and administration for human blood and blood components.

- Gravity blood administration sets
- Blood bag
- Cord blood bag

## Urinary care




Page 43-46

We provide solutions to help improve patient safety and infection prevention in urinary care.

- Urine drainage bag
- Leg bag
- Urine meter
- Irrigation tray
- Irrigation syringe

## Renal therapy



Page 47-48

Innovative, high performance products and services for use in dialysis, provide a basis which targets improved clinical outcomes for the patient, while minimizing costs and risks for the healthcare providers.

- Fistula needle
- Safety fistula needle
- Blood line
- Hemodialysis care kit

## Diabetes care



Page 49-52

Complete range of the diabetes needles and the innovative safety insulin syringe offers comfort and easy insulin therapy Injection, also help eliminate accidental needle stick injury.

- Safety insulin syringe
- Insulin syringe
- Insulin pen needle
- Safety insulin pen needle
- Sharp disposal container for insulin pen needle
- Alcohol swab





## KB MEDICAL(GROUP),INC.

### Our Headquarter

Add.: 9650 Research Drive, Suite A, Irvine, California 92618, USA

Tel: 1-949-333-3785

Email: [info@kbmed.com](mailto:info@kbmed.com)

Website: [www.kbmed.com](http://www.kbmed.com)

### Our production unit

Jiangsu Kangbao Medical Equipment Co. Ltd.

Add.: No. 78, North Suzhong Road, Baoying, Yangzhou city, Jiangsu, China 225800

### Our commercial office

Hangzhou Best Valley Medical Technology Co. Ltd.

Add.: Room 1105, Tower A, Wufutianxinglong Mansion,

No. 358 East Fengqi Road, Jianggan district, Hangzhou, China 310020

Tel: 0086-571-87987813, 87975379, 87983593

Email: [susan@bestvalleymedical.com](mailto:susan@bestvalleymedical.com)

Website: [www.kbmed.com](http://www.kbmed.com)



## KB Medical Group

The Specialist in  
Safety Medical Products  
and Nutrition Care

[www.kbmed.com](http://www.kbmed.com)

PRODUCT CATALOGUE



# KB Medical Group



KB Medical Group was founded in Irvine, California.

A professional company specialized in producing innovative Safety Medical Products to protect health care workers against sharps injury. KB Medical Group is committed to address these issues by providing the most effective safety-engineered therapy devices for optimum safety in the hospitals.

In the past several years, KB Medical Group focus on working on a wide variety of innovative, high-quality enteral feeding products and accessories uniquely designed for delivering nutrition and medication to pediatric and adult patients. This new standard ENFit® connection has been introduced in-line with ISO80369-3 has a unique enteral-specific design that help reduce the risk of enteral tube feeding misconnections and improve the patient safety.


## KB Medical Group in the world

- Our Headquarter: Irvine, California, USA.
- Our production unit: Jiangsu Kangbao Medical Equipment Co. Ltd. is our professional manufacturer in Jiangsu province, China (1 hour away Shanghai) for more than 30 years.
- Our commercial office: Hangzhou Best Valley Medical Technology Co. Ltd. is our Commercial Center in China which is located at Hangzhou, Zhejiang province (40minutes away Shanghai), the center in charge of the business communication with customers, arrange the customers for visiting plan, QA audit, also provide service for the registration documents and documentary works.

## Key facts of our KB Medical Group



## Quality, Certificate, Innovation





## Key facts of our production unit



40 sets

Injection molding machines

16 sets

Auto printing lines for syringes

8 sets

Auto-assembling machines for syringes

8 sets

Blister packing machines

8 sets

Auto-loading packing machines for syringe  
(1/3/5/10/20/60/100ml)

1 set

Auto-assembling machine for infusion set

1 set

Auto-assembling machine for scalp vein set

3 sets

Auto-assembling machines  
for hypodermic needle

3 sets

Modernization testing lab for  
Physical, Chemical and Biological

3 sets

EO Sterilization chamber




▲ Disposablesyringeworkshop



▶ Hypodermicneedleworkshop

▼ Injectionmoldingworkshop



Testing center ▶

