

# Enteral Nutrition



# Enteral Syringe

- FDA cleared for direct oral administration and enteral delivery
- The ENFit® syringe tip is compliant with ISO 80369-3 standards
- ENFit® syringes are often used in pharmacies, NICU/ PICU, adult ICU, alternate care, home care, and physician offices.

## ENFit® Single Use Syringe



### Product Feature :

- ENFit® female.
- DEHP Free
- Latex Free
- Low dose
- Amber color is available



Enteral Nutrition

## ENFit® Home Use O-Ring Syringe

- Up to 84 uses
- O-Ring Plunger
- Low resistance
- Amber color is available

Sizes	Low Dose	Single Use	Home Use
1ml	●	ES1	ESH1
3ml	●	ES3	ESH3
5ml	●	ES5	ESH5
10ml		ES10	ESH10
20ml		ES20	ESH20
35ml		ES35	ESH35
60ml		ES60	ESH60
100ml		ES100	N/A





## Thumb Control Ring Syringe

- Designed for maximum control and one-handed administration
- Re-sealable pole bag designed for safe and convenient storage

## Flat Top Syringe



- Visible dosing marks designed for ease of use

Sizes	Enfit Tip	Thumb Control Ring	Flat Top	Code
60ml		•		URTS60
60ml			•	URFS60
60ml	•	•		URES60
60ml	•		•	ES60
100ml	•		•	ES100
100ml		•		URTS100



## The Enfit® Connector

- The ENFit® connector is compliant with the ISO80369-3.
- Supply you transition adaptor from Legacy connector to Enfit. (ENC3)
- ENFit® syringe cap used as a protective cap to temporarily secure any remaining nutrition between feedings. (ENC6)
- Supply you transition adaptor from Enfit to Legacy connector. (ENC10)



# ENFit® Feeding Set

## Feature :

- ENFit® feeding sets are compliant with the ISO 80369-3 to help reduce the risk of misconnections.
- Latex-free and DEHP-free.
- EVA material available.
- G-Gravity use, P-pump use, ENF01-05 available G type and P type.



**ENF05G**

Universal Set



**ENF01G**

ENPLUS spike with cap



**ENF02P**

Feeding bag  
Sizes: 500/1000/1200/1500ml



**ENF03G**

Feed cap adaptor



**ENF04G**

Spike with feeding connector



Enteral Nutrition



## Oral Dispenser

### Feature :

- Special tip design for wrong route administration prevention.
- O-Ring plunger design is the preferred option for smooth and accurate delivery.
- Amber barrel design to protect light-sensitive drug.




O-Ring Plunger

Type 1



Type 2

Sizes	Type 1	Type 2
1ml	ORS1-01	ORS2-01
3ml	ORS1-03	ORS2-03
5ml	ORS1-05	ORS2-05
10ml	ORS1-10	ORS2-10
20ml	ORS1-20	ORS2-20





## ENFit Stepped

### We're on the path to safer connections

The ENFit connector reduces the risk of tube feeding misconnections, helping enhance patient safety. It has a unique enteral-specific design.



- Has a locking feature that secures the safety connection with administration sets and syringes that have ENFit connectors.
- Aligns with global, industry-wide standards specified by ISO and recommended by the FDA.
- Supply you transition adaptor from Enfit to Legacy connector

ENC10