



Specimen Collection



Blood Collection Needle



Type 1 Butterfly type



Type 2 Pen type



Type 3 Hypodermic needle type



Size		Color	Code									
OD(mm)	Gauge		Safety Type		Safety Type with pre-attached holder			Normal Type		Normal Type with pre-attached holder		
			Type 1	Type 3	Type 1	Type 2	Type 3	Type 1	Type 2	Type 1	Type 2	Type 3
0.40	27G	Grey	SBCN27-1	SBCN27-3	SBCN27-1H	SBCN27-2H	SBCN27-3H	BCN27-1	BCN27-2	BCN27-1H	BCN27-2H	BCN27-3H
0.45	26G	Brown	SBCN26-1	SBCN26-3	SBCN26-1H	SBCN26-2H	SBCN26-3H	BCN26-1	BCN26-2	BCN26-1H	BCN26-2H	BCN26-3H
0.50	25G	Orange	SBCN25-1	SBCN25-3	SBCN25-1H	SBCN25-2H	SBCN25-3H	BCN25-1	BCN25-2	BCN25-1H	BCN25-2H	BCN25-3H
0.55	24G	Purple	SBCN24-1	SBCN24-3	SBCN24-1H	SBCN24-2H	SBCN24-3H	BCN24-1	BCN24-2	BCN24-1H	BCN24-2H	BCN24-3H
0.60	23G	Blue	SBCN23-1	SBCN23-3	SBCN23-1H	SBCN23-2H	SBCN23-3H	BCN23-1	BCN23-2	BCN23-1H	BCN23-2H	BCN23-3H
0.70	22G	Black	SBCN22-1	SBCN22-3	SBCN22-1H	SBCN22-2H	SBCN22-3H	BCN22-1	BCN22-2	BCN22-1H	BCN22-2H	BCN22-3H
0.80	21G	Green	SBCN21-1	SBCN21-3	SBCN21-1H	SBCN21-2H	SBCN21-3H	BCN21-1	BCN21-2	BCN21-1H	BCN21-2H	BCN21-3H
0.90	20G	Yellow	SBCN20-1	SBCN20-3	SBCN20-1H	SBCN20-2H	SBCN20-3H	BCN20-1	BCN20-2	BCN20-1H	BCN20-2H	BCN20-3H
1.10	19G	Cream	SBCN19-1	SBCN19-3	SBCN19-1H	SBCN19-2H	SBCN19-3H	BCN19-1	BCN19-2	BCN19-1H	BCN19-2H	BCN19-3H
1.20	18G	Pink	SBCN18-1	SBCN18-3	SBCN18-1H	SBCN18-2H	SBCN18-3H	BCN18-1	BCN18-2	BCN18-1H	BCN18-2H	BCN18-3H

Blood Collection Tube

We provide a full offering of blood collection tubes in both glass and plastic.

Type	Code
Glass blood collection tube	BCT01
Plastic blood collection tube	BCT02
Micro collection tube	BCT03
Non-vacuum blood collection tube	BCT04
Capillary blood collection tube	BCT05



Category	Item	Additive	Tube Material	Closure Cap					Test Item	
				Type C	Tube Size (mm)	Type H	Tube Size (mm)	Color		
Serum Blood Collection Tube	Plain Tube	Plain	Glass/Plastic		13*75		16*100	Red	Clinical Biochemistry, Immunology and Serology Test	
	Pro-coagulation Tube	Clot Activator	Glass/Plastic		13*100					16*100
	Gel & Clot Activator Tube	Gel & Coagulant	Glass/Plastic		13*75					13*100
Whole Blood Collection Tube	EDTA Tube	Sprayed K2EDTA/ Sprayed K3EDTA	Glass/Plastic		13*75		16*100	Purple	Haematology Test (Blood Routine Examination)	
	ESR Tube	3.8% Sodium Citrate Buffer(0.129mol/L)	Glass/Plastic Glass		13*75					8*120
Plasma Blood Collection Tube	Coagulation Tube	3.2% Sodium Citrate Buffer(0.109mol/L)	Glass/Plastic		13*75		10.25*64	Blue	Coagulation Function Test	
	Heparin Tube	Sodium Heparin/ Lithium Heparin	Glass/Plastic		13*75					13*100
	Gel & Heparin Tube	Gel & Sodium Heparin/ Gel & Lithium Heparin	Glass/Plastic		13*75		16*100	Green	Clinical Chemistry for Emergency Treatment, Blood Rheology Test	
			Glass/Plastic		13*100					16*100
	Glucose Tube	Sodium Fluoride & Sodium Heparin/ Sodium Fluoride & EDTA/ Sodium Fluoride & Potassium Oxalate	Glass/Plastic		13*75		13*100	Grey	Glucose and Lactate Test	
EDTA & Gel Tube	Gel & Sprayed K2EDTA/ Gel & Sprayed K3EDTA	Glass/Plastic	13*75		13*100					16*100



Accessories



Tourniquet

Sizes: 1"x18"x0.025"/1"x18"x0.029"/etc.
Latex free, Blue/Orange/Green/Yellow/Pink



Adhesive bandage

Sizes: 72x19mm/76x25mm/Dia. 25mm/etc.
Sizes can be customized



Alcohol Prep Pads

Sizes: 3x6cm/ 4x8cm



Paper Tape

Width: 1.25cm/2.5cm/5cm/7.5cm/10cm



Self-adhesive bandage

Width: 2.5cm/5cm/7.5cm/10cm/15cm



Cotton ball

Gram weight: 0.3-4 g