



Diabetes Care




Insulin Syringe

Syringes for subcutaneous insulin therapy injection.

Feature:

- Premium needle quality and silicon-coated to ease skin penetration and minimize pain.
- High transparency barrel with clearly readable graduated scale for better dose precision.
- No-dead-space.



Code	Unit	Size	Needle size	Cap color
IS03	U-100	0.3ML	25-33GX5/16"-1/2"	orange
IS05	U-40/ U-100	0.5ML	25-33GX5/16"-1/2"	U-40 red/ U-100 orange
IS10	U-40/ U-100	1ML	25-33GX5/16"-1/2"	U-40 red/ U-100 orange

Needle sizes according to customer's requests



Safety Insulin Syringe

Feature :

- Safety insulin syringes with integrated needle
- The innovative safety insulin syringe offers comfort and easy insulin therapy injection.
- Help eliminate accidental needle stick injury.



Type 1

- Pull back the plunger after injection, the locking ring secures the needle in the barrel of the syringe protecting from needle stick injury.

Type 2

- Pull the safety cartridge to withdraw the needle into cartridge after injection.
- Easy to lock by rotating the safety cartridge.

Code	Size	Needle size
Type1/Type2		
ISS1-03 ISS2-03	0.3ML	29-33G*5/16"-1/2"
ISS1-05 ISS2-05	0.5ML	29-33G*5/16"-1/2"
ISS1-10 ISS2-10	1ML	28-33G*5/16"-1/2"





Insulin Pen Needle

Feature:

- Extremely fine three-facet grinding minimizes tissue trauma
- Special silicone finish for smooth gliding and virtually pain-free penetration of skin and tissue
- Extra-thin-wall technology allows for a larger inner diameter and an optimized Insulin flow; less force is required and injection time is reduced
- Compatible with all major pen devices available on the market



Safety Insulin Pen Needle

Feature:

- Automatically achieve the function of protecting needle tip after use.
- Easy self-disable set
- Easy to see the droplets
- Compatible with all pen injectors



Code	Gauze	size and length (mm)
Safety/Non-safety		
IPNS33 IPN33 new	33G	0.20X4-12mm
IPNS32 IPN32	32G	0.23X4-12mm
IPNS31 IPN31	31G	0.25X4-12mm
IPNS30 IPN30	30G	0.3X4-12mm
IPNS29 IPN29	29G	0.33x4-12mm

Product Features



Safety design and easy for using

The strength design in needle set make it more convenient and safe to use

High compatibility and cover almost all branches
The product has high compatibility. The standard design of internal threaded ensure it can stably match all branches insulin pen.


Best penetration make injection more comfortable
Insulin pen needle has a series of advanced and mature needle tip technique and silication technique. The sharp needle tip make puncture more easily and let patients feel more comfortable during injection.



Diabetes Care



33G, a shorter and finer insulin pen needle!!



Don't need to pinch the skin when use 4mm needle

Contrast of 8mm, 4mm needle used in few subcutaneous tissue area

High Matching Degree



NovoPen® 4



NovoPen® 3



NovoPen® Junior



NovoRapid® Flexpen®



BD: Puregon Pen®



Sanofi: Aventis Lantus



Sanofi: OptiPen®



Sanofi: TactiPen®



Owen Mumford: AutoPen®



Owen Mumford: AutoPen®2



Lily: HumaPen®



Lily: HumaPen®



Wanbang: Wangbang Pen 2



Tong hua dong bao: Doing Bao Pen®



Tong hua dong bao: Gansulin Pen®

- ▶ The screw thread design of insulin pen needle meets the standard regulated by ISO11608-2 which is compatible with other standard ones domestically produced or imported.



Safety Lancet

It is used in blood collection sampling and help eliminate accidental needlestick injuries.

Feature:

- Needle is fully shielded before and after use
- Cannot be reused
- Ultra sharp tri-bevel needle, high speed penetration
- Multiple depth selection with color coding



Type3



Type1



Type4



Type5

Code					Size	Length
TYPE 1	TYPE 2	TYPE 3	TYPE 4	TYPE 5		
SBL01-30	SBL02-30	SBL03-30	SBL04-30	SBL05-30	30G	1.8mm
SBL01-28	SBL02-28	SBL03-28	SBL04-28	SBL05-28	28G	1.8mm
SBL01-26	SBL02-26	SBL03-26	SBL04-26	SBL05-26	26G	1.8mm
SBL01-23	SBL02-23	SBL03-23	SBL04-23	SBL05-23	23G	1.8mm
SBL01-21	SBL02-21	SBL03-21	SBL04-21	SBL05-21	21G	1.8/2.2/2.8mm
SBL01-18	SBL02-18	SBL03-18	SBL04-18	SBL05-18	18G	1.8mm

Blood Glucose Monitoring Systems

Feature :

- Fast, Safe and Convenient.
- Small blood sample requirement reduces pain & sensitivity for easy diabetes monitoring.
- Self-testing on the volume of the blood sample
- Coding/no-coding optional



Small blood sample
Less pain



Alarm



5 seconds
Fast testing



Auto-off



Strip ejector



High accuracy
ISO 15197:2013
certified



No coding



Test strip shelf
life 2 yearsTest



360 test memories



BG205

Feature :

- 1ul small blood sample
- 5seconds measuring time
- 360 reading memories
- Extra large LCD
- Backlight optional
- Strip ejector
- 7,14,28 day average
- Alarm clock reminder
- Switchable dual scale
- Auto Power off
- Deluxe carry case
- Coding/no-coding optional



Diabetes Care



BG202

Feature :

- 1ul small blood sample
- 5 seconds measuring time
- 360 reading memories
- 7,14,28 day average
- Alarm clock reminder
- Switchable dual scale
- Auto Power off
- Deluxe carry case
- Coding/ no-coding optional



BG203

Feature :

- 1ul small blood sample
- 5 seconds measuring time
- 360 reading memories
- Extra large LCD
- Backlight optional
- Strip ejector
- 7,14,28 day average
- Alarm clock reminder
- Switchable dual scale
- Auto Power off
- Deluxe carry case
- Coding/ no-coding optional

BG208

Feature :

- 1ul small blood sample
- 5 seconds measuring time
- Strip ejector
- 360 reading memories
- 7,14,28 day average
- Alarm clock reminder
- Switchable dual scale
- Auto Power off
- Deluxe carry case
- Coding/nocoding optional



Accessories



Sharp disposal container

0.2L/1Quart/1L/2Quart/ 2L



Alcohol swab

3x6cm/ 4x8cm



Lancet

23g /26g /28g /30g /32g /33g



Lancing devices

The adjustable tip offer 5 depth levels



Safety Blood Lancet

18g/21g/23g /26g /28g /30g



Test strips

Width 5.5mm/6.5mm

