

# Transfusion therapy



## Gravity blood administration sets



Filtration for various blood IV administration needs.

### Features

- ▶ Sharp spike is easy to puncture.
- ▶ Plastic spike/metallic spike available.
- ▶ Vented with hydrophobic bacterial filter.
- ▶ The ISO standard spike can fit with any blood bags.
- ▶ Soft and transparent chamber.
- ▶ Extra soft & kink resistant tubes.
- ▶ Acute roller clamp.
- ▶ With 200micron filtration.
- ▶ Latex free.
- ▶ Straight and Y type are available.
- ▶ DEHP free tube is available.



BTS001

● Big chamber with air inlet



BTS002

● Big chamber without air inlet



BTS003

● Blow chamber

# Transfusion therapy



## Blood bag

It is used to collect and preserve whole human blood and blood components in hospital or blood bank.



### Features

- ▶ Our blood bag could not only guarantee high quality of the blood components but also to provide maximum safety during blood donation and sample collection.
- ▶ The safety needle shield enhances safety against the accidental needle stick injury that can transmit HIV, Hepatitis B and C and other blood borne diseases.
- ▶ Closed system & reusable blood sampling with vacuum tube.
- ▶ Label material with pre-cutting line and could not be tear integrated.



**BB01/Single**  
Single blood bag is used for collection of whole human blood. It consists of one bag containing anticoagulant solution CPDA.

**BB02/Double**  
Double bag system consists of one primary bag containing anticoagulant solution CPD/CPDA and one empty satellite bag, for separation of half plasma.

**BB03/Triple**  
Triple bag system consists of one primary bag containing anticoagulant solution CPD/CPDA, one satellite bag containing SAGM preservation solution, and one empty satellite bag. This additive can maintain red blood cell viability, with over 80% of survival rate after 5 weeks of storage.

**BB04/Quadruple**  
Quadruple bag system consists of one primary bag containing anticoagulant solution CPD/CPDA, one satellite bag containing SAGM preservation solution, and two empty satellite bags. The additive SAGM can maintain red blood cell viability, with over 80% of survival rate after 5 weeks of storage.



BB03/Triple

# Transfusion therapy



## Cord blood bag

It is designed to collect the umbilical cord blood (also known as placental blood) after the birth of the baby, within the shortest time and minimal risk of contamination. The clinical use of the cord blood is allergenic transplantation in patients suffering from severe blood diseases such as leukemia, lymphoma, etc.



CBB01



Leukocyte filter



Sampling device & safety needle protector

### Optional

- ▶ Top and top system.
- ▶ Top and bottom system.
- ▶ Empty transfer bag.
- ▶ PALL leukocyte filtration system.
- ▶ Blood sampling with vacuum tube.
- ▶ Customization bags.

# CATALOGUE CATEGORIES



## Infusion therapy



Page 1-9

- Wide range of a safety way for injection of vein for medication infusion.
- Safety scalp vein set
  - Scalp vein set
  - Infusion set
  - IV burette
  - IV cannula
  - Safety IV cannula
  - IV cannula dressing
  - Sterile IV start kit
  - Extension tube
  - Three way stopcock

## Syringe and needle



Page 10-21

A complete range of syringes for injection or harvesting of biological material with safety protection.

- Disposable syringe with needle
- Disposable syringe with safety needle
- Safety syringe
- Safety needle
- Disposable hypodermic needle
- Auto-Disable syringe for fixed-dose immunization
- Pre-filled flush syringe
- Dental needle
- Irrigation needle
- Suture removal Set
- Veterinary needle
- Hard pack syringe

## Hazardous drug safety



Page 22-28

Various drug preparation needles for safety drug admixture and disposal of used sharps.

- Aspiration needle with safety shield
- Blunt needle with/without filter
- Curved needle
- Side hole needle
- Needle free system
- Safety box for the disposal of used sharps
- Sharp container
- Endoscope transport pad
- Ultimate sharps safety stations (needle counter)
- Safety medical products for NEEDLE STICK PREVENTION

## Safety medical products



Page 28

Wide range of safety medical products to protect health care workers against sharps injuries.

- Safety scalp vein set
- Safety IV cannula
- Safety needle
- Safety syringe
- Auto-Disable syringe
- Safety blood collection needle
- Safety blood lancet
- Safety insulin pen needle
- Safety insulin syringe
- Safety fistula needle
- Safety drug preparation needle
- Ultimate sharps safety stations
- Needle free system
- Safety box for the disposal of used sharps
- Sharp container

## Nutrition therapy



Page 29-32

Various components for nutrition therapy which is compliant with the ISO80369-3 to help reduce the risk of misconnections.

- ENFit@syringe
- The ENFit@ connector
- ENFit@ Feeding set
- Oral syringe
- Irrigation syringe

## Specimen collection



Page 33-36

Our safety blood collection system help eliminate accidental needle stick injuries when collection blood sample.

- Blood collection needle
- Safety blood collection needle
- Safety lancet
- Blood collection tube

## Anesthesia & pain management



Page 37-39

Wide range of products for subdural, lower thorax and lumbar spinal puncture.

- Spinal needle
- Anesthesia kit
- High pressure manifold
- NRFit™ syringe

## Transfusion therapy



Page 40-42

Complete range of safety way to collect, preserve, filtration and administration for human blood and blood components.

- Gravity blood administration sets
- Blood bag
- Cord blood bag

## Urinary care



Page 43-46

We provide solutions to help improve patient safety and infection prevention in urinary care.

- Urine drainage bag
- Leg bag
- Urine meter
- Irrigation tray
- Irrigation syringe

## Renal therapy



Page 47-48

Innovative, high performance products and services for use in dialysis, provide a basis which targets improved clinical outcomes for the patient, while minimizing costs and risks for the healthcare providers.

- Fistula needle
- Safety fistula needle
- Blood line
- Hemodialysis care kit

## Diabetes care



Page 49-52

Complete range of the diabetes needles and the innovative safety insulin syringe offers comfort and easy insulin therapy Injection, also help eliminate accidental needle stick injury.

- Safety insulin syringe
- Insulin syringe
- Insulin pen needle
- Safety insulin pen needle
- Sharp disposal container for insulin pen needle
- Alcohol swab





## KB MEDICAL(GROUP),INC.

### Our Headquarter

Add.: 9650 Research Drive, Suite A, Irvine, California 92618, USA

Tel: 1-949-333-3785

Email: [info@kbmed.com](mailto:info@kbmed.com)

Website: [www.kbmed.com](http://www.kbmed.com)

### Our production unit

Jiangsu Kangbao Medical Equipment Co. Ltd.

Add.: No. 78, North Suzhong Road, Baoying, Yangzhou city, Jiangsu, China 225800

### Our commercial office

Hangzhou Best Valley Medical Technology Co. Ltd.

Add.: Room 1105, Tower A, Wufutianxinglong Mansion,

No. 358 East Fengqi Road, Jianggan district, Hangzhou, China 310020

Tel: 0086-571-87987813, 87975379, 87983593

Email: [susan@bestvalleymedical.com](mailto:susan@bestvalleymedical.com)

Website: [www.kbmed.com](http://www.kbmed.com)



## KB Medical Group

The Specialist in  
Safety Medical Products  
and Nutrition Care

[www.kbmed.com](http://www.kbmed.com)

PRODUCT CATALOGUE





## Key facts of our production unit



40 sets

Injection molding machines

16 sets

Auto printing lines for syringes

8 sets

Auto-assembling machines for syringes

8 sets

Blister packing machines

8 sets

Auto-loading packing machines for syringe  
(1/3/5/10/20/60/100ml)

1 set

Auto-assembling machine for infusion set

1 set

Auto-assembling machine for scalp vein set

3 sets

Auto-assembling machines  
for hypodermic needle

3 sets

Modernization testing lab for  
Physical, Chemical and Biological

3 sets

EO Sterilization chamber




▲ Disposablesyringeworkshop



▶ Hypodermicneedleworkshop



▼ Injectionmoldingworkshop



Testing center ▶





# KB Medical Group



KB Medical Group was founded in Irvine, California.

A professional company specialized in producing innovative Safety Medical Products to protect health care workers against sharps injury. KB Medical Group is committed to address these issues by providing the most effective safety-engineered therapy devices for optimum safety in the hospitals.

In the past several years, KB Medical Group focus on working on a wide variety of innovative, high-quality enteral feeding products and accessories uniquely designed for delivering nutrition and medication to pediatric and adult patients. This new standard ENFit® connection has been introduced in-line with ISO80369-3 has a unique enteral-specific design that help reduce the risk of enteral tube feeding misconnections and improve the patient safety.

## KB Medical Group in the world

- Our Headquarter: Irvine, California, USA.
- Our production unit: Jiangsu Kangbao Medical Equipment Co. Ltd. is our professional manufacturer in Jiangsu province, China (1 hour away Shanghai) for more than 30 years.
- Our commercial office: Hangzhou Best Valley Medical Technology Co. Ltd. is our Commercial Center in China which is located at Hangzhou, Zhejiang province (40minutes away Shanghai), the center in charge of the business communication with customers, arrange the customers for visiting plan, QA audit, also provide service for the registration documents and documentary works.

## Key facts of our KB Medical Group



## Quality, Certificate, Innovation

