

# Enteral Syringe

- FDA cleared for direct oral administration and enteral delivery
- The ENFit® syringe tip is compliant with ISO 80369-3 standards
- ENFit® syringes are often used in pharmacies, NICU/ PICU, adult ICU, alternate care, home care, and physician offices.

## ENFit® Single Use Syringe



### Product Feature :

- ENFit® female.
- DEHP Free
- Latex Free
- Low dose
- Amber color is available



Enteral Nutrition

## ENFit® Home Use O-Ring Syringe

- Up to 84 uses
- O-Ring Plunger
- Low resistance
- Amber color is available

Sizes	Low Dose	Single Use	Home Use
1ml	●	ES1	ESH1
3ml	●	ES3	ESH3
5ml	●	ES5	ESH5
10ml		ES10	ESH10
20ml		ES20	ESH20
35ml		ES35	ESH35
60ml		ES60	ESH60
100ml		ES100	N/A



# ENFit® Feeding Set

## Feature:

- ENFit® feeding sets are compliant with the ISO 80369-3 to help reduce the risk of misconnections.
- Latex-free and DEHP-free.
- EVA material available.
- G-Gravity use, P-pump use, ENF01-05 available
- G type and P type.



**ENF05G**

Universal Set



**ENF01G**

ENPLUS spike with cap



**ENF02P**

Feeding bag  
Sizes: 500/1000/1200/1500ml



**ENF03G**

Feed cap adaptor



**ENF04G**

Spike with feeding connector



Enteral Nutrition



## ENFit Stepped

### We're on the path to safer connections

The ENFit connector reduces the risk of tube feeding misconnections, helping enhance patient safety. It has a unique enteral-specific design.



- Has a locking feature that secures the safety connection with administration sets and syringes that have ENFit connectors.
- Aligns with global, industry-wide standards specified by ISO and recommended by the FDA.
- Supply you transition adaptor from Enfit to Legacy connector

ENC10