



*Since 1986*



# Blood Collection Tube



## **KB Medical Group Inc.**

We Focus on a Strong Customer  
Orientation and Innovative  
Top Quality Product Solutions.

5325 Alton Parkway Ste C 732, Irvine, CA 92604








T: +1 (949) 333-3785    E: info@kbmed.com

[www.kbmed.com](http://www.kbmed.com)







## Disposable Vacuum Blood Tube

Category	Item	Additive	Type Cap	Volume	Tube Size (mm)	Test items	Code					
							Glass Tube	PET Tube				
Plasma Blood Collection Tube				10ml	16*100		BCT01EK3100	BCT02EK3100				
				8ml			BCT01EK3080H	BCT02EK3080H				
				9ml			BCT01EK3090H	BCT02EK3090H				
				10ml			BCT01EK3100H	BCT02EK3100H				
	ESR Tube	3.8% Sodium Citrate		1.6ml	13*75	Erythrocyte Sedimentation Rate Test	BCT01ESR016	BCT02ESR016				
				2.4ml			BCT01ESR024	BCT02ESR024				
				3.2ml			BCT01ESR032	BCT02ESR032				
				3.6ml			BCT01ESR036	BCT02ESR036				
				4.0ml	13*100		BCT01ESR040	BCT02ESR040				
				4.8ml			BCT01ESR048	BCT02ESR048				
				5.6ml	8*120		BCT01ESR056	BCT02ESR056				
				1.28ml			BCT01ESR0128H					
				1.6ml			BCT01ESR016H					
				1.6ml			BCT01ESRN016H					
				Coagulation Tube	3.2% Sodium Citrate			1.8ml	13*75	Coagulation Function Test	BCT01SC018	BCT02SC018
								2.7ml			BCT01SC027	BCT02SC027
								3.6ml			BCT01SC036	BCT02SC036
								4.5ml			BCT01SC045	BCT02SC045
	5.4ml	13*100	BCT01SC054			BCT02SC054						
	6.3ml		BCT01SC063			BCT02SC063						
	Heparin Tube	Sodium Heparin				1ml		13*75	Clinical Chemistry for Emergency Treatment, Blood Rheology Test		BCT01SH010	BCT02SH010
						2ml					BCT01SH020	BCT02SH020
						3ml					BCT01SH030	BCT02SH030
						4ml					BCT01SH040	BCT02SH040
				5ml	13*100	BCT01SH050	BCT02SH050					
				6ml		BCT01SH060	BCT02SH060					
				7ml	16*100	BCT01SH070	BCT02SH070					
				8ml		BCT01SH080	BCT02SH080					
9ml				BCT01SH090		BCT02SH090						
10ml				BCT01SH100		BCT02SH100						
Lithium Heparin			1ml	13*75	BCT01LH010	BCT02LH010						
			2ml		BCT01LH020	BCT02LH020						
			3ml		BCT01LH030	BCT02LH030						
			4ml		BCT01LH040	BCT02LH040						
			5ml	13*100	BCT01LH050	BCT02LH050						
			6ml		BCT01LH060	BCT02LH060						
			7ml	16*100	BCT01LH070	BCT02LH070						
			8ml		BCT01LH080	BCT02LH080						
			9ml		BCT01LH090	BCT02LH090						
			10ml		BCT01LH100	BCT02LH100						
Gel & Heparin Tube	Gel & Sodium Heparin		1ml	13*75	BCT01SHG010	BCT02SHG010						
			2ml		BCT01SHG020	BCT02SHG020						
			3ml		BCT01SHG030	BCT02SHG030						
			4ml		BCT01SHG040	BCT02SHG040						
			5ml	13*100	BCT01SHG050	BCT02SHG050						
			6ml		BCT01SHG060	BCT02SHG060						
			7ml	16*100	BCT01SHG070	BCT02SHG070						
			8ml		BCT01SHG080	BCT02SHG080						
			8.5ml		BCT01SHG085	BCT02SHG085						
			1ml		13*75	BCT01LHG010	BCT02LHG010					
	2ml	BCT01LHG020	BCT02LHG020									
	3ml	BCT01LHG030	BCT02LHG030									
	4ml	BCT01LHG040	BCT02LHG040									
	Gel & Lithium Heparin		5ml	13*100	BCT01LHG050	BCT02LHG050						
			6ml		BCT01LHG060	BCT02LHG060						
			7ml	16*100	BCT01LHG070	BCT02LHG070						
			8ml		BCT01LHG080	BCT02LHG080						
			8.5ml		BCT01LHG085	BCT02LHG085						
							BCT01SFH010	BCT02SFH010				
					BCT01SFH020	BCT02SFH020						

## Disposable Vacuum Blood Tube

Category	Item	Additive	Type Cap	Volume	Tube Size (mm)	Test items	Code			
							Glass Tube	PET Tube		
Glucose Tube	Sodium Fluoride & Sodium Heparin			3ml	13*75	Glucose and Lactate Test	BCT01SFH030	BCT02SFH030		
				4ml			BCT01SFH040	BCT02SFH040		
				5ml			BCT01SFH050	BCT02SFH050		
				6ml	13*100		BCT01SFH060	BCT02SFH060		
				7ml			BCT01SFH070	BCT02SFH070		
				1ml			BCT01SNA010	BCT02SNA010		
				2ml	13*75		BCT01SNA020	BCT02SNA020		
				3ml			BCT01SNA030	BCT02SNA030		
				4ml			BCT01SNA040	BCT02SNA040		
				5ml	13*100		BCT01SNA050	BCT02SNA050		
				6ml			BCT01SNA060	BCT02SNA060		
				7ml			BCT01SNA070	BCT02SNA070		
				1ml	13*75			BCT01SP010	BCT02SP010	
				2ml				BCT01SP020	BCT02SP020	
				3ml				BCT01SP030	BCT02SP030	
	4ml	13*100	BCT01SP040	BCT02SP040						
	5ml		BCT01SP050	BCT02SP050						
	6ml		BCT01SP060	BCT02SP060						
	7ml	BCT01SP070	BCT02SP070							
	EDTA & Gel Tube	Gel & K2 EDTA			1ml	13*75	Molecular Biology Test	BCT01EK2G010	BCT02EK2G010	
					2ml			BCT01EK2G020	BCT02EK2G020	
					3ml			BCT01EK2G030	BCT02EK2G030	
					4ml	13*100		BCT01EK2G040	BCT02EK2G040	
					5ml			BCT01EK2G050	BCT02EK2G050	
					6ml			BCT01EK2G060	BCT02EK2G060	
					7ml	16*100		BCT01EK2G070	BCT02EK2G070	
					8ml			BCT01EK2G080	BCT02EK2G080	
					8.5ml			BCT01EK2G085	BCT02EK2G085	
					1ml	13*75			BCT01EK3G010	BCT02EK3G010
					2ml				BCT01EK3G020	BCT02EK3G020
					3ml				BCT01EK3G030	BCT02EK3G030
					4ml	13*100		BCT01EK3G040	BCT02EK3G040	
					5ml			BCT01EK3G050	BCT02EK3G050	
					6ml			BCT01EK3G060	BCT02EK3G060	
		7ml	16*100	BCT01EK3G070	BCT02EK3G070					
		8ml		BCT01EK3G080	BCT02EK3G080					
8.5ml		BCT01EK3G085		BCT02EK3G085						

## Micro Blood Collection Tube

Code	Additive	Tube Material	Type Color	Volume	Tube Size (mm)
BCT03A010	Plain	PP	Red	1ml	10x45
BCT03A005				0.5ml	
BCT03A0025				0.25ml	
BCT03CA010	Clot Activator	PP	Red	1ml	10x45
BCT03CA005				0.5ml	
BCT03CA0025				0.25ml	
BCT03EK2010	EDTA K2	PP	Purple	1ml	10x45
BCT03EK2005				0.5ml	
BCT03EK20025				0.25ml	
BCT03EK3010	EDTA K3	PP	Purple	1ml	10x45
BCT03EK3005				0.5ml	
BCT03EK30025				0.25ml	
BCT03LH010	Heparin Lithium	PP	Green	1ml	10x45
BCT03LH005				0.5ml	
BCT03LH0025				0.25ml	
BCT03SH010	Heparin Sodium	PP	Green	1ml	10x45
BCT03SH005				0.5ml	
BCT03SH0025				0.25ml	
BCT03SFH010	Heparin Sodium & Sodium Fluoride	PP	Grey	1ml	10x45
BCT03SFH005				0.5ml	
BCT03SFH0025				0.25ml	
BCT03SGC008	Gel & Clot Activator	PP	Yellow	0.8ml	10x45
BCT03SGC005				0.5ml	
BCT03SGC0025				0.25ml	

