


Anesthesia & pain management



High pressure manifold



Maximize fluid flow control and fluid and drug delivery efficiency. Use two and three port normally closed valve manifolds for maximum fluid and drug delivery efficiency in the Operating Room. Confidently maintain drug dosing integrity when administering multiple drugs with back check valves in each port.



Features

- ▶ No leakage under 500psi.
- ▶ Smooth rotation, easy use.
- ▶ Different color for choice.
- ▶ DEHP-free, Latex-free.

Code	Spec.
HPM2	Dual Valve
HPM3	Triple Valve
HPM4	Four Valve

NRFit™ syringe

It is used for regional anesthesia, the new patient safety initiative neuraxial products with NRFit™ connectors according to American Society of Anesthesiologists (ASA) meeting, Oct. 22-26 in Chicago 2016. All healthcare facilities are expected to begin converting to meet the ISO standard requiring NRFit connectors for neuraxial use. Healthcare facilities has already taken legislative action by mandating that convert to these new standardized epidural products by Jan. 1, 2017.



Features

- ▶ NRFit™ connectors design-- ISO80369-6 compliant.
- ▶ Neuraxial applicant connector.
- ▶ Small bore connectors that are intended to prevent misconnections to improve patient safety.

Code	Size	Type
NRS3	3ml	NRFit™ connector
NRS5	5ml	NRFit™ connector
NRS10	10ml	NRFit™ connector
NRS20	20ml	NRFit™ connector

Spinal needle



Subdural, lower thorax and lumbar spinal puncture.

Features

- ▶ Full sizes of anesthesia needles.
- ▶ Identified as spinal needle bevel quincke tip, pencil point tip and epidural needle.
- ▶ 304 stainless steel needle tubes and stylets.
- ▶ International colour coding.
- ▶ With mounted introducer needle.
- ▶ Luer lock hub.



Quincke tip

Quincke tip	Introducer	Length	Code
18G×3-1/2"	None	38-220mm	SNQ18
19G×3-1/2"	None	38-220mm	SNQ19
20G×3-1/2"	None	38-220mm	SNQ20
21G×3-1/2"	None	38-220mm	SNQ21
22G×3-1/2"	None	38-220mm	SNQ22
23G×3-1/2"	19G×1-1/2"	38-220mm	SNQ23
24G×3-1/2"	20G×1-1/2"	38-220mm	SNQ24
25G×3-1/2"	20G×1-1/2"	38-220mm	SNQ25
26G×3-1/2"	20G×1-1/2"	38-220mm	SNQ26
27G×3-1/2"	22G×1-1/2"	38-220mm	SNQ27

- ▶ Trauma-free QUINCKE short tip to ensure immediate injection point closure and reduced CSF loss.
- ▶ Differentiated angle coupling of needle and obturator tip to avoid production and movement of tissue fragments, improving operator sensitivity.



Pencil point

Pencil point	Introducer	Length	Code
22g×3-1/2"	None	90-220mm	SNP22
23G×3-1/2"	19G×1-1/2"	90-220mm	SNP23
24G×3-1/2"	20G×1-1/2"	90-220mm	SNP24
25G×3-1/2"	20G×1-1/2"	90-220mm	SNP25
26G×3-1/2"	20G×1-1/2"	90-220mm	SNP26
27G×3-1/2"	22G×1-1/2"	90-220mm	SNP27

- ▶ With no sharp edges separates fibers without slicing.
- ▶ The lateral hole shape and its position relative to the tip favour rapid CSF reflux and effective analgesic administration.
- ▶ Reduce the risk of saddle block.

Anesthesia kit



Single-use anesthesia set is clinically used for epidural or spinal anesthesia during puncture and drug injection. This product can only be for single-use, in order to avoid cross-infection. It contains the unique components needed to perform anesthesia operation.



AK01

Components

- ▶ Epidural needle, Size: 15G,16G,17G,18G. Length: 80mm, 90mm.
- ▶ Spinal needle, Type: Quincke tip & pencil point, Size: 18G,19G,20G,21G,22G,23G,24G,25G,26G,27G. Length: 38mm–150mm.
- ▶ Epidural catheter, Size: OD: ϕ 1.0, OD: ϕ 0.8, Length: 850mm.
- ▶ Epidural filter: 0.22 μ m with hydrophobic membrane luer lock.
- ▶ LOR syringe: 5ml, 7ml, 10ml.
- ▶ Catheter adaptor: TBA, secure the catheter completely.

Optional

- ▶ Epidural kit, spinal kit and combined kit available.
- ▶ Catheter with /without metal wire available.
- ▶ Precision filter 0.2-5 μ m available.
- ▶ Double layer UV filter 0.2-5 μ m are available.
- ▶ Any customization components.

CATALOGUE CATEGORIES



Infusion therapy



Page 1-9

- Wide range of a safety way for injection of vein for medication infusion.
- Safety scalp vein set
 - Scalp vein set
 - Infusion set
 - IV burette
 - IV cannula
 - Safety IV cannula
 - IV cannula dressing
 - Sterile IV start kit
 - Extension tube
 - Three way stopcock

Syringe and needle



Page 10-21

A complete range of syringes for injection or harvesting of biological material with safety protection.

- Disposable syringe with needle
- Disposable syringe with safety needle
- Safety syringe
- Safety needle
- Disposable hypodermic needle
- Auto-Disable syringe for fixed-dose immunization
- Pre-filled flush syringe
- Dental needle
- Irrigation needle
- Suture removal Set
- Veterinary needle
- Hard pack syringe

Hazardous drug safety



Page 22-28

Various drug preparation needles for safety drug admixture and disposal of used sharps.

- Aspiration needle with safety shield
- Blunt needle with/without filter
- Curved needle
- Side hole needle
- Needle free system
- Safety box for the disposal of used sharps
- Sharp container
- Endoscope transport pad
- Ultimate sharps safety stations (needle counter)
- Safety medical products for NEEDLE STICK PREVENTION

Safety medical products



Page 28

Wide range of safety medical products to protect health care workers against sharps injuries.

- Safety scalp vein set
- Safety IV cannula
- Safety needle
- Safety syringe
- Auto-Disable syringe
- Safety blood collection needle
- Safety blood lancet
- Safety insulin pen needle
- Safety insulin syringe
- Safety fistula needle
- Safety drug preparation needle
- Ultimate sharps safety stations
- Needle free system
- Safety box for the disposal of used sharps
- Sharp container

Nutrition therapy



Page 29-32

Various components for nutrition therapy which is compliant with the ISO80369-3 to help reduce the risk of misconnections.

- ENFit@syringe
- The ENFit@ connector
- ENFit@ Feeding set
- Oral syringe
- Irrigation syringe

Specimen collection



Page 33-36

Our safety blood collection system help eliminate accidental needle stick injuries when collection blood sample.

- Blood collection needle
- Safety blood collection needle
- Safety lancet
- Blood collection tube

Anesthesia & pain management



Page 37-39

Wide range of products for subdural, lower thorax and lumbar spinal puncture.

- Spinal needle
- Anesthesia kit
- High pressure manifold
- NRFit™ syringe

Transfusion therapy



Page 40-42

Complete range of safety way to collect, preserve, filtration and administration for human blood and blood components.

- Gravity blood administration sets
- Blood bag
- Cord blood bag

Urinary care



Page 43-46

We provide solutions to help improve patient safety and infection prevention in urinary care.

- Urine drainage bag
- Leg bag
- Urine meter
- Irrigation tray
- Irrigation syringe

Renal therapy



Page 47-48

Innovative, high performance products and services for use in dialysis, provide a basis which targets improved clinical outcomes for the patient, while minimizing costs and risks for the healthcare providers.

- Fistula needle
- Safety fistula needle
- Blood line
- Hemodialysis care kit

Diabetes care



Page 49-52

Complete range of the diabetes needles and the innovative safety insulin syringe offers comfort and easy insulin therapy Injection, also help eliminate accidental needle stick injury.

- Safety insulin syringe
- Insulin syringe
- Insulin pen needle
- Safety insulin pen needle
- Sharp disposal container for insulin pen needle
- Alcohol swab



KB MEDICAL(GROUP),INC.

Our Headquarter

Add.: 9650 Research Drive, Suite A, Irvine, California 92618, USA

Tel: 1-949-333-3785

Email: info@kbmed.com

Website: www.kbmed.com

Our production unit

Jiangsu Kangbao Medical Equipment Co. Ltd.

Add.: No. 78, North Suzhong Road, Baoying, Yangzhou city, Jiangsu, China 225800

Our commercial office

Hangzhou Best Valley Medical Technology Co. Ltd.

Add.: Room 1105, Tower A, Wufutianxinglong Mansion,

No. 358 East Fengqi Road, Jianggan district, Hangzhou, China 310020

Tel: 0086-571-87987813, 87975379, 87983593

Email: susan@bestvalleymedical.com

Website: www.kbmed.com



KB Medical Group

The Specialist in
Safety Medical Products
and Nutrition Care

www.kbmed.com

PRODUCT CATALOGUE



KB Medical Group



KB Medical Group was founded in Irvine, California.

A professional company specialized in producing innovative Safety Medical Products to protect health care workers against sharps injury. KB Medical Group is committed to address these issues by providing the most effective safety-engineered therapy devices for optimum safety in the hospitals.

In the past several years, KB Medical Group focus on working on a wide variety of innovative, high-quality enteral feeding products and accessories uniquely designed for delivering nutrition and medication to pediatric and adult patients. This new standard ENFit® connection has been introduced in-line with ISO80369-3 has a unique enteral-specific design that help reduce the risk of enteral tube feeding misconnections and improve the patient safety.


KB Medical Group in the world

- Our Headquarter: Irvine, California, USA.
- Our production unit: Jiangsu Kangbao Medical Equipment Co. Ltd. is our professional manufacturer in Jiangsu province, China (1 hour away Shanghai) for more than 30 years.
- Our commercial office: Hangzhou Best Valley Medical Technology Co. Ltd. is our Commercial Center in China which is located at Hangzhou, Zhejiang province (40minutes away Shanghai), the center in charge of the business communication with customers, arrange the customers for visiting plan, QA audit, also provide service for the registration documents and documentary works.

Key facts of our KB Medical Group



Quality, Certificate, Innovation



Key facts of our production unit



40 sets

Injection molding machines

16 sets

Auto printing lines for syringes

8 sets

Auto-assembling machines for syringes

8 sets

Blister packing machines

8 sets

Auto-loading packing machines for syringe
(1/3/5/10/20/60/100ml)

1 set

Auto-assembling machine for infusion set

1 set

Auto-assembling machine for scalp vein set

3 sets

Auto-assembling machines
for hypodermic needle

3 sets

Modernization testing lab for
Physical, Chemical and Biological

3 sets

EO Sterilization chamber



▲ Disposablesyringeworkshop



▶ Hypodermicneedleworkshop

▼ Injectionmoldingworkshop



Testing center ▶

