

The premier manufacturer of ENFit® products for optimal enteral nutrition.

Why choose KB Medical for your nutrition feeding systems? As we are the leading global manufacturer of safety-engineered therapy devices for optimal hospital and patient safety. We go above and beyond to give you exceptional quality, service and value.

Safety Support

As a medical safety products manufacturer focused on protecting health care workers against sharps injury, we expanded our safety to Enteral Feeding Products.

We manufacture ISO80369-3 compliant enteral feeding products and accessories uniquely designed using the new standard ENFit® connection system for effective delivery of nutrition and medication to pediatric and adult without misconnection.

Manufacturer-Direct Pricing

- Unlike other providers, we control the entire supply chain.
- Eliminate third-party costs.

Ensure the Top Quality

- Growing network of ISO-certified and quality-compliant manufacturing centers.
- Rigorous quality control and testing.
- Reliable and trusted products for clinicians and caregivers to optimize clinical efficiency.

Customized production satisfy different customer's requirement

We could provide OEM brand products and innovation technology design to fulfill different customer's requirements.

- Customized packing
- Customized design
- Customized orders

Hands-on Service and Ease of conversion

Having a readily available representative and support network with the expertise to problem solve with you is priceless.

- Proactive attention of a dedicated Representative focused on your success.
- Our experts in Irvine, California USA will work with you to ensure the conversion is smooth and timely.



ENTERAL FEEDING

GRAVITY SET



KB Medical (Group) Inc.

9650 Research Drive, Suite A, Irvine California 92618 USA

www.kbmednutrition.com | www.kbmed.com

T: 1-949-333-3785 F: 1-949-409-0206 Email: info@kbmed.com

Safely and Completely Reliable

When your nutrition delivery system is user friendly and dependable, everyone benefits.

- Nurses are more efficient.
- Patients receive better care.
- Families are more confident in your facility.

Our KBMED™ Enteral Feeding Gravity Sets have features that make operation and formula administration easy and safe.

Key features

- Enfit and Enfit transition connectors are available
- Regular and large bores are available
- Bags are available in multiple size 500ml/1000ml/1200ml

Multiple Components



The Enfit connector is compliant with the ISO80369-3 to help reduce the risk of misconnections.



Large bore flow regulator flow rate 3times faster than regular ones.



Large bore tube
OD: 6.0MM
ID: 4.1MM

Regular bore tube
OD: 4.1MM
ID: 3.0MM

Main materials

- Medical grade PVC tubing.
- DEHP free and latex free.

We're on the path to safer connections

The ENFit connector reduces the risk of tube feeding misconnections, helping enhance patient safety. It has a unique enteral-specific design that:

- Has a locking feature that secures the safety connection with administration sets and syringes that have ENFit connectors.
- Aligns with global, industry-wide standards specified by ISO and recommended by the FDA.
- Supply you transition adaptor from Enfit to Legacy connector



We're available with Large-bore Enteral Feeding Gravity Sets

- KBMED™ Enteral Feeding Gravity sets are available with large bore.
- Gastrostomies use a large-bore tube, which allows more viscous feedings and decrease risk of tube occlusion.

ITEM CODE	FEEDING BAG SIZE			DRIP CHAMBER	FLOW REGULATOR		TUBE DIAMETER		PATIENT LABEL	ENFIT CONNECTOR	ENFIT TRANSITION CONNECTOR
	500ML	1000ML	1200ML		REGULAR	LARGE	4.1*3MM	6.0*4.1MM			
ENF02G-1000R-T		•		•	•		•		•	•	•
ENF02G-1000B-T		•		•		•		•	•	•	•
ENF02G-1200R-T			•	•	•		•		•	•	•
ENF02G-1200B-T			•	•		•		•	•	•	•
ENF02G-500R-T	•			•	•		•		•	•	•
ENF02G-1000R		•		•	•		•		•	•	
ENF02G-1000B		•		•		•		•	•	•	
ENF02G-1200R			•	•	•		•		•	•	
ENF02G-1200B			•	•		•		•	•	•	
ENF02G-500R	•			•	•		•		•	•	