

# Urinary Care



# Urine Drainage Bag

Urinary drainage bags are intended to collect urine from sheaths or catheters. It provides solutions to help improve patient safety and infection prevention in urinary care.

## Features:

- Anti-reflux drip chamber on the bag helps prevent urine reflux.
- With airvent and non-return valve
- Varies choice of components: universal hanger, bed sheet clamp, dog house, needless samples port.
- High grade PVC.
- Latex free.
- Tyvek pouch for O.R. pack or PE packing.



Urinary Care



### UDBA20

Apple shape 2000ml  
4000ml is available.



### UDBD20

Diamond shape 2000ml



## Urine Meter

### Features :

- Designed for easy bag suspension with a reinforced double hanger and built -in anti-back guide.
- Convenient needleless sample port is accessible with either needleless or blunt cannula, and offers added protection for patient and caregiver.
- Adjustable bed sheet clamp secures tubing.
- Include additional rope hanger.



Code	Spec.
UM01	400ml meter with 2600ml bag, 3 chambers.
UM02	400ml meter with 2000ml bag, 2 chambers.
UM03	500ml meter with 1500ml bag, 4 chambers.
UM04	200ml meter with 2000ml bag, 1 chamber.
UM05	500ml meter with 2000ml bag, 1 chamber.
Remark:200ml-500ml capacity meter is available	

### UDBRT20

Rectangle shape with T tap outlet 2000ml



### UDBRP20

Rectangle shape With pull outlet 2000ml



### UDBR20

Rectangle shape without outlet 2000ml

## Leg Bag

- Rectangle or arc style for customer option.
- It is secured with two comfortable straps in latex free, they have a twist drainage outlet or flip flow valve for easy emptying.
- T-tap outlet or twist outlet for one hand easy emptying.
- With white fabric for comfortable skin touch.
- Tube length 30-90cm for customer option.

Code	Spec.
LB500	500ml
LB600	600ml
LB750	750ml
LB900	900ml
LB1000	1000ml



Urinary Care

## Irrigation Tray

Irrigation tray is designed to improve simplicity and efficiency of care.

### Components

Qty	Description
1ea	1000ml collection tray
1ea	500 ml graduated container
1ea	60ml bulb or thumb control ring syringe
1ea	Protector cap
2ea	Alcohol prep pad
1ea	Moisture-proof under pad 45x45cm
1pr	Vinyl gloves



URIT1000